



Introducing the respirator that exceeds N95 performance



95%
Particulate
protection
efficiency

New KN95 Protective Mask

Approved under the FDA EUA

Specifications	N95 NIOSH	KN95 GB2626-2006
• Filter Performance (must be greater \geq X% efficient)	$\geq 95\%$	$\geq 95\%$ (Actual 99.21% -99.4%)
• Test Agent	NaCl	NaCl
• Flow rate	85 L/min	85 L/min
• Total inward leakage (TIL) tested on human subjects each performing exercises	N/A	$\leq 8\%$ leakage (arithmetic mean)
• Inhalation Resistance (maximum pressure drop)	≤ 343 Pa	≤ 350 Pa
• Flow rate	85 L/min	85 L/min
• Exhalation Resistance (maximum pressure drop)	≤ 245 Pa	≤ 250 Pa
• Flow rate	85 L/min	85 L/min